

**Saint Luke's Hospital
Kansas City, MO 64111**

**Thromboelastograph (TEG) Request Form
Call Blood Bank (22412) before sending sample.
Send to pneumatic tube station #12.**

Date/Time of collection (use time on computer) _____ **A#** _____

Is patient on Heparin and/or Protamine ?

Yes (Blood Bank run a clear and heparinase cup)

No (Blood Bank run a clear cup only)

Specimen Collection

- Draw blood into blue top citrated tube to the fill line.
- Do not open the cap to fill the tube.
- Mix by gentle inversion five times.
- A separate specimen **MUST** be drawn for a COAG screen.

Line Draw (10 mL waste drawn prior to sample draw)

Venipuncture (3 mL waste drawn prior to sample draw)